



# Carbon Nanotubes For Next Generation X-Ray Systems

Moritz Beckmann, XinRay Systems LLC



National Nanotechnology Initiative at Ten:  
**NANOTECHNOLOGY**  
**INNOVATION** SUMMIT

**DECEMBER 8-10, 2010**  
GAYLORD CONVENTION CENTER  
WASHINGTON, DC



# XinRay Systems – the company





7020 Kit Creek Rd, Suite 210  
 Research Triangle Park  
 NC 27709  
<http://www.xinraysystems.com>

- Founded in Aug 2007 as JV between Xintek and Siemens Medical Solutions
  - 11 highly skilled employees
  - Located in Research Triangle Park, NC
- 
- Experts and equipment from Siemens and Xintek
  - Development contracts in security and medical



# Carbon Nanotubes as Field Emitters





CNTs work like inverse lightning rods

# XinRay Technology



## Conventional technology





Carbon Nanotubes (CNTs) are “cold emitters”

- Faster than thermionic emitters
- Cold → less energy required
- Can be placed close together

- Individual emitters of multibeam array
- can be activated in rapid sequence
  - Mimic movement
  - Multiplexing possible

## XinRay x-ray tube technology



→ Enables 3D imaging systems without moving parts



# Advantages – Very fast switching – Less Dose





## Preclinical applications:

- Periodic mechanical shutter
- Out of sync data is disregarded

## XinRay solution:

- Extra fast switching of x-ray
- Electronically controlled exposures
- Non – periodic synchronization
- Less dose




# Advantages – No mechanical movement – fast imaging





## Digital Breast Tomosynthesis:

- Improved detection
- But too long exam time  
→ patient discomfort and motion blur



## Stationary Digital Breast Tomosynthesis:

- Faster imaging
- No motion blur






# Advantages – No mechanical movement – fast imaging



## Image Guidance in Radiation Therapy:

- Irradiation of tumor from different directions
- Several sittings to treat cancer
- Planning CT, then positioning of patient relative to markers





# Advantages – No mechanical movement – fast imaging






## Image Guidance in Radiation Therapy:

- Irradiation of tumor from different directions
- Several sittings to treat cancer
- Planning CT, then positioning of patient relative to mark

## Tomosynthesis (limited angle tomography):

- Acquisition of image in a few seconds
- Imaging right before and even during treatment (interfractional imaging)
- High resolution of tumor position in plane perpendicular to treatment beam.








# Shortcomings of conventional CT

For 40 years in 3D x-ray imaging a single source (+ detector) is moved around object to image from different directions

low resolution



insufficient detection



manual check



>85.5 million bags in 2006  
(16%)

**Cost of False Alarms - \$ 1 billion / year**



## Disadvantages – conventional technology

- Limitation in speed → limited throughput
- Mechanical wear and tear → high downtime and maintenance costs
- Single x-ray source → motion blur & lower resolution → high false alarm rate
- System geometry dominated by circular movement → large space requirements

**Slow - High false alarm - Bulky - High maintenance**



## **X**Advantages – XinRay technology

- **No** Limitation in speed → **4X throughput increase**
- **No** Mechanical wear and tear → **lower downtime and maintenance costs**
- **Better resolution → better detection** → **lower false alarm rate**
- **Geometry optimized for application** → **smaller / better detection**

**This technology “will revolutionize the CT industry”.**

Dr. Zarur, Science Advisor, Department of Homeland Security (DHS), in Nature, July 2009





With Carbon Nanotubes XinRay offers  
**Faster – Smaller - Better imaging – More robust**  
Solutions and also enables  
**New Generation of Imaging Systems**

**Thank You !**

